

Fungal Morphogenesis

Fungal structures such as hyphae, spores, fruiting bodies, and mycelial networks are critical to a wide range of biological and ecological processes. These include nutrient recycling, symbiosis, pathogenicity, and environmental resilience. Their analysis is essential for understanding fungal biology, species identification, ecological interactions, and the diagnosis of fungal infections in both clinical and environmental settings.

There is strong commercial and medical interest in fungal imaging, particularly for managing agricultural diseases, diagnosing human and animal infections, and developing biotechnological and pharmaceutical products such as antibiotics, immunosuppressants, and industrial enzymes. In pathogenic fungi like *Candida albicans*, morphological transitions are directly linked to virulence, drug resistance, and host invasion, making detailed structural analysis essential for both research and therapeutic monitoring.



Conventional light microscopy can reveal general fungal morphology but often lacks the resolution to distinguish fine structural details, especially in dense, overlapping, or heterogeneous samples. Features such as septation, early hyphal emergence, and internal compartments may go undetected. While higher-resolution methods can overcome these limitations, they typically involve complex workflows that can obscure native structural relationships. These limitations highlight the need for imaging methods that combine high resolution, minimal preparation, and preservation of native structure.

Morphological Diversity in Fixed *Candida albicans* Populations


We have observed that fixed populations of *Candida albicans* display a mixture of yeast-like cells, pseudo-hyphal chains, and early hyphal forms within the same sample. Rounded cells are frequently seen clustered together, likely representing recently budded individuals or spore-like structures.

Chains of elongated cells with visible constrictions suggest pseudo-hyphal growth—a transitional morphology often associated with increased invasive capacity. Some cells exhibit polarised extensions at one end, characteristic of the onset of hyphal development. These elongated forms frequently display distinct internal contrast, which may correspond to vacuoles, nuclei, or intracellular storage granules, providing insight into the cells' metabolic and developmental state.

HalolImages reveal structural features in Fixed *Candida albicans* Populations:



HalolImage reveals multiple stages of fungal development within the same population. The zoomed-in section on the right highlights several elongated cells arranged in short chains, characteristic of pseudo-hyphal growth – a transitional form between budding and true filamentous growth. (Scale bar: 20 μ m).



HaloImage shows diverse morphogenetic stages, with the inset highlighting a rounded, yeast-like cell and an elongated form indicative of early filamentous growth.

(Scale bar: 20 μ m).

Additional zoomed-in sections from the left image reveal additional structures including budding cells, pseudo-hyphal chains, and early hyphal extensions.

Together, these features highlight the dynamic structural variation within *C. albicans* populations and illustrate transitions between yeast, filamentous, and reproductive forms under fixed conditions. By capturing this diversity with clarity and internal detail, HaloMicroscopy supports a deeper understanding of fungal morphogenesis and pathogenesis.

HaloMicroscopy also extends to live-cell imaging. Its ability to track morphological transitions—such as budding, hyphal extension, and spore germination—in real time within hydrated environments enables applications in antifungal drug screening, biofilm development studies, and dynamic, time-lapse analysis of fungal behaviour.